

## **CARING FOR YOUR BITE APPLIANCE**

Proper care of your appliance is critical to maintain its proper fit, function and to ensure its longevity.

- Rinse your appliance after every use with warm water and clean with a soft toothbrush. Do not use toothpaste to clean your appliance. Use only a non abrasive cleaner such as a denture approved cleaning solution or mild dish detergent
- Store your appliance DRY in its supplied case
- Keep your storage case clean
- Do not clean your appliance with boiling water as it will result in distortion. Use only warm water that is not too hot for your hands
- Any issues with your appliance, please call the office.

We are here to help if you need us! Please contact our office at 610-667-2770 or [info@drweinstock.com](mailto:info@drweinstock.com) if you have any questions or issues with these instructions.